





**Learn more about column packing!**  
Watch our instructional column packing videos available on YouTube™



**Make sure that you are using the right chromatography medium (resin):**

Download the free Purify App from GE Healthcare's Life Sciences

[www.gelifesciences.com/purify](http://www.gelifesciences.com/purify)



**Focus on your research by letting GE pack the column for you!**

With more than 50 years of experience in chromatography and column packing, GE offers:

- Pre-packed columns
- Custom column packing

[www.gelifesciences.com/custom-column-packing](http://www.gelifesciences.com/custom-column-packing)

Product name	Max. operating pressure			Column inner diameter (mm)	Approx. bed volume (ml)	Bed height (mm)
	bar	psi	MPa			
Tricorn 5/20 Column	100	1450	10	5	0.0-0.5	0-26
Tricorn 5/50 Column	100	1450	10	5	0.2-1.1	8-56
Tricorn 5/100 Column	100	1450	10	5	1.2-2.1	58-106
Tricorn 5/150 Column	100	1450	10	5	2.1-3.1	108-156
Tricorn 5/200 Column	100	1450	10	5	3.1-4.1	158-206
Tricorn 10/20 Column	50	725	5	10	0.0-2.1	0-26
Tricorn 10/50 Column	50	725	5	10	0.0-4.4	0-56
Tricorn 10/100 Column	50	725	5	10	3.6-8.4	46-106
Tricorn 10/150 Column	50	725	5	10	7.6-12.3	96-156
Tricorn 10/200 Column	50	725	5	10	11.5-16.2	146-206
Tricorn 10/300 Column	50	725	5	10	19.4-24.1	246-306
Tricorn 10/600 Column	50	725	5	10	42.9-47.6	546-606
HiScale 16/20 Column	20	290	2	16	0-40	0-200
HiScale 16/40 Column	20	290	2	16	16-80	80-400
HiScale 26/20 Column	20	290	2	26	0-106	0-200
HiScale 26/40 Column	20	290	2	26	69-212	130-400
HiScale 50/20 Column	20	290	2	50	0-393	0-200
HiScale 50/40 Column	20	290	2	50	274-785	140-400
XK 16/20 Column	5	72.5	0.5	16	0-31	0-155
XK 16/40 Column	5	72.5	0.5	16	16-70	80-350
XK 16/70 Column	5	72.5	0.5	16	76-130	380-650
XK 16/100 Column	5	72.5	0.5	16	136-190	680-950
XK 26/20 Column	5	72.5	0.5	26	0-66	0-125
XK 26/40 Column	5	72.5	0.5	26	45-186	85-350
XK 26/70 Column	5	72.5	0.5	26	204-344	385-650
XK 26/100 Column	5	72.5	0.5	26	363-504	685-950
XK 50/20 Column	3	43.5	0.3	50	0-274	0-140
XK 50/30 Column	3	43.5	0.3	50	0-559	0-280
XK 50/60 Column	3	43.5	0.3	50	500-1088	260-560
XK 50/100 Column	3	43.5	0.3	50	1274-1862	650-950
C 10/10 Column	1	14.5	0.1	10	0-7	0-80
C 10/20 Column	1	14.5	0.1	10	0-14	0-180
C 10/40 Column	1	14.5	0.1	10	16-30	200-380
C 16/20 Column	1	14.5	0.1	16	0-36	0-180
C 16/40 Column	1	14.5	0.1	16	8-77	40-380
C 16/70 Column	1	14.5	0.1	16	69-137	348-680
C 16/100 Column	1	14.5	0.1	16	129-197	640-980
C 26/40 Column	1	14.5	0.1	26	21-202	40-380
C 26/70 Column	1	14.5	0.1	26	181-361	348-680
C 26/100 Column	1	14.5	0.1	26	340-520	640-980
Empty Disposable PD-10 Column	-	-	-	-	≤ 13	-